NADAL[®] Dengue NS1 Ag + IgG/IgM Test (test cassette)

REF 532019N-10

de	Gebrauchsanweisung	2	Symbols	39
en	Instructions for use	7	Our Teams	40
fr	Instructions d'utilisation	12		
es	Instrucciones de uso	17		
it	Istruzioni per l'uso	22		
pl	Sposób użycia	27		
pt	Instruções de Utilização	32		

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14. Referenzen

- Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. Rev. Infect. Dis. 1984; 6:251-264.
- Halstead SB, Pathogenesis of dengue: challenges to molecular biology. Science 1988; 239:476-481.
- Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). Japanese J. Trop. Med. Hygiene 1994; 22: 9-12.
- Lam SK. Dengue haemorrhagic fever. Rev. Med. Micro. 1995; 6:39-48.
- Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
- Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. Clin Diagn Lab Immunol. 2003 Jul; 10(4): 725-8.
- Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. Wien Med Wochenschr (in German). 1997; 147(19-20): 463-4.
- Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus. Microbiol Immunol. 1994; 38(12): 951-5.

Rev. 0, 2020-10-19 OM

1. Intended Use

The NADAL® Dengue NS1 Ag + IgG/IgM Test is a lateral flow chromatographic immunoassay for the qualitative detection of dengue virus NS1 antigen as well as anti-dengue virus IgG and IgM in human whole blood, serum or plasma specimens (see section 12 'Limitations'). This test is intended for use as an aid in the diagnosis of primary and secondary dengue infections. The test procedure is not automated and requires no special training or qualification. The NADAL® Dengue NS1 Ag + IgG/IgM Test is designed for professional use only.

2. Introduction and Clinical Significance

Dengue virus is a flavivirus transmitted primarily by tiger mosquitos *Aedes aegypti* and *Aedes albopictus*. It is widely distributed throughout the tropical and subtropical areas of the world and causes up to 100 million infections annually. Classic dengue infection is characterised by a sudden onset of fever, intense headache, myalgia, arthralgia and a rash.

Primary dengue infection causes IgM to increase to a detectable level approximately 5 days after the onset of fever. IgM generally persists for 30 to 90 days. Most dengue patients in endemic regions have secondary infections, resulting in high levels of specific IgG prior to or simultaneous with IgM response. Therefore, the detection of specific anti-dengue IgM and IgG can also help to distinguish between primary and secondary infections.

NS1 is one of 7 dengue virus non-structural proteins which are thought to be involved in viral replication. NS1 exists as a monomer in its immature form but is rapidly processed in the endoplasmic reticulum to form a stable dimer. A small amount of NS1 remains associated with intracellular organelles, where it is thought to be involved in viral replication. The rest of NS1 is found to be either associated with the plasma membrane or secreted as a soluble hexadimer. NS1 is essential for viral viability but its precise biological function is unknown. Antibodies produced in response to NS1 in viral infections can cross-react with cell surface antigens on epithelial cells and platelets. This may have implications in the development of dengue haemorrhagic fever.

3. Test Principle

The NADAL® Dengue NS1 Ag + IgG/IgM Test is a lateral flow chromatographic immunoassay for the qualitative detection of dengue virus NS1 antigens as well as anti-dengue virus IgG and IgM in human whole blood, serum or plasma specimens.

The NADAL® Dengue NS1 Ag + IgG/IgM Test contains two internal test strips (left side of the test cassette: Dengue NS1 Ag Test; right side of the test cassette: Dengue IgG/IgM Test).

Anti-dengue virus NS1 antibodies are immobilised in the test line region 'T' of the membrane. Anti-human IgM and anti-human IgG are immobilised in the test line regions 'IgM' and 'IgG', respectively, of the membrane.

During the test, dengue virus NS1 antigens, if present in the specimen, bind to anti-dengue virus NS1 antibodies (active ingredient of the detection reaction: 0.6 - 0.8 µg/test), which are conjugated to coloured particles and precoated onto the sample pad of the test cassette. If specific anti-dengue virus IgM or anti-dengue virus IgG are present in the specimen, they bind to the dengue virus antigens (active ingredient of the

detection reaction: 0.2 - 0.4 µg/test) conjugated to coloured particles and precoated onto the sample pad of the test cassette. The mixture then migrates along the membrane by capillary action and interacts with the reagents on the membrane. The complexes are then captured by anti-dengue virus NS1 antibodies (active ingredient of the detection reaction: 1.0 - 2.5 µg/test) as well as anti-human IgM (active ingredient of the detection reaction: 0.5 - 1.0 µg/test) and anti-human IgG (active ingredient of the detection reaction: 1.5 - 3.0 µg/test) in the test line regions 'T' as well as 'IgM' and 'IgG', respectively. Excess coloured particles are captured in each control line region 'C'. The presence of a coloured line in the test line region 'T' and/or 'IgM' and/or 'IgG' indicates a positive result. The absence of a coloured line in the test line region 'T' and/or 'IgM' and/or 'IgG' indicates a negative result.

The formation of a coloured line in the control line region 'C' (left side of the test cassette: Dengue NS1 Ag Test) and the colour change from a blue line to a red one in the control line region 'C' (right side of the test cassette: Dengue IgG/IgM Test) serve as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 10 NADAL® Dengue NS1 Ag + IgG/IgM test cassettes*, incl. disposable pipettes (25 µL) for Dengue NS1 Ag
- 10 disposable pipettes (5 µL) for Dengue IgG/IgM
- 1 buffer (3 mL)**
- 1 package insert

*containing the preservative sodium azide: <0.1%

**Phosphate buffer containing the following preservatives: sodium azide: <0.1%

No hazard labelling is required according to Regulation (EC) № 1272/2008 CLP. Concentrations are below exemption threshold.

5. Additional Materials Required

- Specimen collection containers (appropriate for specimen material to be tested)
- Centrifuge (for serum or plasma specimens only)
- Alcohol pads
- Lancets (for fingerstick whole blood specimens only)
- Timer

6. Storage and Stability

Test kits should be stored at 2-30°C until the indicated expiry date. Test cassettes are stable until the expiry date printed on the foil pouches. Test cassettes must remain in the sealed foil pouches until use. Do not freeze the test kit. Do not use tests beyond the expiry date indicated on the packaging. Care should be taken to protect test kit components from contamination. Do not use test kit components if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to inaccurate results.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure prior to testing.

- Do not use the test beyond the expiration date indicated on the packaging.
- Do not use test kit components if the primary packaging is damaged.
- Tests are for single use only.
- Do not add specimens to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new specimen collection container for each specimen obtained.
- Do not substitute or mix components from different test kits.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Temperature can adversely affect test results.
- Used testing materials should be disposed of according to local regulations.

8. Specimen Collection and Preparation

The NADAL® Dengue NS1 Ag + IgG/IgM Test can be performed using whole blood (from venipuncture or fingerstick), serum or plasma.

To collect fingerstick whole blood specimens:

- Wash the patient's hand with soap and warm water or clean it with an alcohol pad. Allow it to dry.
- Massage the hand, without touching the puncture site, by rubbing along the hand towards the fingertip of the middle or ring finger.
- Puncture the skin with a sterile lancet. Wipe away the first drop of blood.
- Gently rub the hand from the wrist to the palm, and then to the finger to form a rounded drop of blood over the puncture site.

Fingerstick whole blood should be tested immediately.

Venipuncture whole blood specimens

Containers containing anticoagulants, such as sodium EDTA, sodium citrate or heparin sodium should be used for the preparation of venous whole blood or plasma specimens.

Testing should be performed immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods of time.

If the test is to be run within 3 days of specimen collection, whole blood collected by venipuncture should be stored at 2-8°C.

Do not freeze whole blood specimens.

Serum and plasma specimens

Separate serum or plasma from blood as soon as possible to avoid haemolysis. Use only clear, non-haemolysed specimens.

Testing should be performed immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods of time. Serum and plasma specimens can be stored at 2-8°C for up to 7 days. For long-term storage, specimens should be kept at -20°C.

Bring specimens to room temperature prior to testing. Frozen specimens should be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.

If specimens are to be shipped, they should be packed in compliance with all applicable regulations for the transportation of etiologic agents.

Icteric, lipemic, haemolysed, heat-treated and contaminated specimens may lead to inaccurate test results.


9. Test Procedure

Bring tests, specimens, buffer and/or controls to room temperature (15-30°C) prior to testing.

1. Remove the test cassette from the foil pouch and use it as soon as possible. The best results will be obtained if the test is performed immediately after opening the foil pouch. Label the test cassette with the patient or control identification.
2. Place the test cassette on a clean and level surface.
3. **a) For serum or plasma specimens:**


Dengue virus NS1 Ag

Holding the 25 µL-pipette vertically, draw up a serum or plasma specimen and add 3 drops (approximately 75 µL) of the specimen to the specimen well (S1) of the test cassette. **Avoid air bubbles forming.**

Anti-dengue virus IgG/IgM


- Holding the 5 µL-pipette vertically, draw a serum or plasma specimen up to the first widening (approximately 5 µL) and add it to the specimen well (S2) of the test cassette.
- Holding the buffer bottle vertically, add 2 drops of buffer to the specimen well (S2). **Avoid air bubbles forming.**

b) For venepuncture/fingerstick whole blood specimens:

Dengue virus NS1 Ag

- Holding the 25 µL-pipette vertically, draw up a whole blood specimen and add 2 drops (approximately 50 µL) of the specimen to the specimen well (S1) of the test cassette.


- Holding the buffer bottle vertically, add 1 drop of buffer to the specimen well (S1). **Avoid air bubbles forming.**

Anti-dengue virus IgG/IgM

- Holding the 5 µL-pipette vertically, draw a whole blood specimen up to above the first widening and add 1 drop (approximately 10 µL) of the specimen to the specimen well (S2) of the test cassette.
- Holding the buffer bottle vertically, add 2 drops of buffer to the specimen well (S2). **Avoid air bubbles forming.**

Alternatively, a micropipette may be used.

4. Start the timer.
5. Wait for the red line(s) to appear. Read the test result after exactly 15 minutes. Do not interpret the result after more than 30 minutes.

Negative for anti-dengue virus IgG and IgM



The blue line in the control line region 'C' turns red. No lines develop in the test line region 'IgM' and 'IgG'.

Invalid for dengue virus NS1 Ag

The control line 'C' fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded.

Invalid for anti-dengue virus IgG and IgM

Results from any test where the blue line in the control line region 'C' remains completely or partially blue at the specified reading time must be discarded.

Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line failure.

11. Quality Control

An internal procedural control is included in the test cassette: A coloured line appearing in the control line region 'C' (left side of the test cassette: Dengue NS1 Ag Test) as well as a blue line turning red (right side of the test cassette: Dengue IgG/IgM Test) are considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Good laboratory practice (GLP) recommends the use of external control materials to ensure proper test kit performance.

12. Limitations

- The NADAL® Dengue NS1 Ag + IgG/IgM Test is for professional *in-vitro* diagnostic use only. It should be used for the qualitative detection of dengue virus NS1 antigens as well as anti-dengue virus IgG and IgM in human whole blood, serum or plasma specimens only. Neither the quantitative value nor the rate of increase in the concentration of dengue virus NS1 antigens or anti-dengue virus antibodies can be determined with this qualitative test.
- The NADAL® Dengue NS1 Ag + IgG/IgM Test only detects the presence of dengue virus NS1 antigens and anti-dengue virus antibodies in specimens and should not be used as the sole criterion for a diagnosis of dengue.
- Shortly after the onset of fever, anti-dengue virus IgM concentrations may be below detectable levels. For primary infection, an IgM antibody-capture enzyme-linked immunosorbent assay (MAC-ELISA) showed that 80% of the tested dengue patients exhibited detectable levels of IgM by the 5th day after infection and 99% of the patients tested

10. Result Interpretation


Positive for dengue virus NS1 Ag

A coloured line develops in the control line region 'C' and another coloured line develops in the test line region 'T'.


Positive for anti-dengue virus IgM

The blue line in the control line region 'C' turns red. Another red line develops in the test line region 'IgM'.


Positive for anti-dengue virus IgG

The blue line in the control line region 'C' turns red. Another red line develops in the test line region 'IgG'.

Positive for anti-dengue virus IgG and IgM

The blue line in the control line region 'C' turns red. A red line develops in the test line region 'IgM' and another one develops in the test line region 'IgG'.

Note: The colour intensity in the test line regions 'T', 'IgG' or 'IgM' may vary depending on the concentration of dengue virus NS1 antigens or anti-dengue virus antibodies in the specimen. Therefore, any shade of colour in the test line region 'T', 'IgG' or 'IgM' should be considered positive. Note that this is a qualitative test only and it cannot determine the analyte concentration in the specimen.

Negative for dengue virus NS1 Ag

A coloured line develops in the control line region 'C'. No line develops in the test line region 'T'.

IgM positive by the tenth day. Therefore, it is recommended that patients be tested within this period of time. As for secondary infection, a low molar fraction of anti-dengue virus IgM and a high molar fraction of anti-dengue virus IgG that are broadly reactive to flaviviruses characterise the antibodies. The test line 'IgM' may be faint, and a cross-reaction in the test line region 'IgG' may occur.

- Serological cross-reactivity across the flavivirus group (St. Louis encephalitis, West Nile virus, Japanese encephalitis and yellow fever viruses) is possible. Positive results should be confirmed using other methods.
- The continued presence or absence of antibodies cannot be used to determine the success or failure of therapy.
- Results from immunosuppressed patients should be interpreted with caution.
- As with all diagnostic tests, all results should be interpreted in conjunction with other clinical information available to the physician.
- If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of a recent dengue infection.
- A negative test result for dengue virus NS1 antigens can occur if the quantity of dengue virus NS1 antigens present in the specimen is below the detection limit of the test or if the detected dengue virus NS1 antigens were not present during the stage of disease in which the specimen was collected.
- The presence of detectable dengue virus NS1 Ag may mean a positive result for early dengue infection.
- Dengue virus NS1 antigen is expected to be detectable 1 day after the onset of fever and may persist for up to 9 days in both primary and secondary dengue infections. But if anti-dengue virus NS1 antibodies are produced, the detection of dengue virus NS1 antigens is inhibited.

13. Performance Characteristics

Clinical performance

Diagnostic sensitivity and specificity

The NADAL® Dengue NS1 Ag + IgG/IgM Test was evaluated with clinical specimens whose status was confirmed using ELISA for the right side of the test cassette (Dengue IgG/IgM Test) and using viral culture for the left side of the test cassette (Dengue NS1 Ag Test).

The results are presented in the following tables:

NADAL® Dengue NS1 Ag + IgG/IgM Test (IgM) vs. ELISA

NADAL® Dengue NS1 Ag + IgG/IgM Test (IgM)	Dengue IgM ELISA			
	Positive	Negative	Total	
	Positive	59	0	59
	Negative	1	140	141
Total	60	140	200	

Diagnostic sensitivity: 98.3% (91.1% - 99.7%)*

Diagnostic specificity: >99.9% (97.3% - 100.0%)*

Overall agreement: 99.5% (97.2% - 99.9%)*

*95% confidence interval

NADAL® Dengue NS1 Ag + IgG/IgM Test (IgG) vs. ELISA

NADAL® Dengue NS1 Ag + IgG/IgM Test (IgG)	Dengue IgG ELISA			
	Positive	Negative	Total	
	Positive	89	0	89
	Negative	1	150	151
Total	90	150	240	

Diagnostic sensitivity: 98.9% (94.0% - 99.8%)*

Diagnostic specificity: >99.9% (97.5% - 100.0%)*

Overall agreement: 99.6% (97.7% - 99.9%)*

*95% confidence interval

NADAL® Dengue NS1 Ag + IgG/IgM Test (NS1 Ag) vs. viral culture

NADAL® Dengue NS1 Ag + IgG/IgM Test (NS1 Ag)	Viral culture			
	Positive	Negative	Total	
	Positive	104	3	107
	Negative	8	186	194
Total	112	189	301	

Diagnostic sensitivity: 92.9% (86.5% - 96.3%)*

Diagnostic specificity: 98.4% (95.4% - 99.5%)*

Overall agreement: 96.3% (93.6% - 97.9%)*

*95% confidence interval

Interfering substances

Negative, low positive anti-dengue virus IgG and IgM plasma specimens as well as low positive dengue virus NS1 Ag plasma specimens spiked with the following interfering substances showed no interference with the NADAL® Dengue NS1 Ag + IgG/IgM Test.

Acetylsalicylic acid	4.34 mmol/L
Ascorbic acid	0.2 mg/mL
Albumin	60 g/L
Bilirubin	342 µmol/L
Caffeine	308 µmol/L
Ethanol	86.8 mmol/L
Haemoglobin	2 g/L
HAMA	780.73 IU/mL
Ibuprofen	2425 µmol/L
Metronidazole	701 µmol/L
Paracetamol	199 µmol/L
Quinine	148 µmol/L
Rheumatoid factors	1035 IU/mL
Rifampicin	78.1 µmol/L
Triglycerides	37 mmol/L
γ globulin	60 g/L

Cross-reactivity

Anti-HIV, HBsAg, anti-HCV, anti-HAV IgM, anti-rubella virus IgM, anti-CMV IgM, anti-HSV-1 IgM, anti-HSV-2 IgM, typhoid fever, yellow fever, Zika virus, chikungunya virus, anti-Toxoplasma gondii IgM, syphilis and tuberculosis positive specimens were tested using the NADAL® Dengue NS1 Ag + IgG/IgM Test. No cross-reactivity with the specimens was observed when tested using the NADAL® Dengue NS1 Ag + IgG/IgM Test.

Precision

Repeatability and reproducibility

Precision was established by testing 10 replicates of 7 specimens: a negative, a low and a high positive NS1 Ag

specimen as well as low and high positive anti-Dengue virus IgG and IgM specimens.

Repeatability was established by testing 7 specimens in 3 assays: a negative, a low and a high positive NS1 Ag specimen as well as low and high positive anti-dengue virus IgG and IgM specimens. Testing was performed by 3 operators using 3 independent NADAL® Dengue NS1 Ag + IgG/IgM test lots at 3 different sites on 3 separate days.

The NADAL® Dengue NS1 Ag + IgG/IgM Test demonstrated acceptable repeatability and reproducibility. The negative and positive values were correctly identified >99% of the time.

14. References

1. Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. *Rev. Infect. Dis.* 1984; 6:251-264.
2. Halstead SB, Pathogenesis of dengue: challenges to molecular biology. *Science* 1988; 239:476-481.
3. Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). *Japanese J. Trop. Med. Hygiene* 1994; 22: 9-12.
4. Lam SK. Dengue haemorrhagic fever. *Rev. Med. Micro.* 1995; 6:39-48.
5. Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
6. Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. *Clin Diagn Lab Immunol.* 2003 Jul; 10(4): 725-8.
7. Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. *Wien Med Wochenschr* (in German). 1997; 147(19-20): 463-4.
8. Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus. *Microbiol Immunol.* 1994; 38(12): 951-5.

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14. Bibliographie

- Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. *Rev. Infect. Dis.* 1984; 6:251-264.
- Halstead SB, Pathogenesis of dengue: challenges to molecular biology. *Science* 1988; 239:476-481.
- Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). *Japanese J. Trop. Med. Hygiene* 1994; 22: 9-12.
- Lam SK. Dengue haemorrhagic fever. *Rev. Med. Micro.* 1995; 6:39-48.

- Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
- Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. *Clin Diagn Lab Immunol.* 2003 Jul; 10(4): 725-8.
- Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. *Wien Med Wochenschr (in German).* 1997; 147(19-20): 463-4.
- Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus. *Microbiol Immunol.* 1994; 38(12): 951-5.

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14. Referencias

1. Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. *Rev. Infect. Dis.* 1984; 6:251-264.
2. Halstead SB, Pathogenesis of dengue: challenges to molecular biology. *Science* 1988; 239:476-481.
3. Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). *Japanese J. Trop. Med. Hygiene* 1994; 22: 9-12.
4. Lam SK. Dengue haemorrhagic fever. *Rev. Med. Micro.* 1995; 6:39-48.
5. Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.

6. Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. *Clin Diagn Lab Immunol.* 2003 Jul; 10(4): 725-8.
7. Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. *Wien Med Wochenschr (in German).* 1997; 147(19-20): 463-4.
8. Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus. *Microbiol Immunol.* 1994; 38(12): 951-5.

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14. Bibliografia

- Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. Rev. Infect. Dis. 1984; 6:251-264.
- Halstead SB, Pathogenesis of dengue: challenges to molecular biology. Science 1988; 239:476-481.
- Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). *Japanese J. Trop. Med. Hygiene* 1994; 22: 9-12.
- Lam SK. Dengue haemorrhagic fever. *Rev. Med. Micro.* 1995; 6:39-48.
- Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
- Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. *Clin Diagn Lab Immunol.* 2003 Jul; 10(4): 725-8.
- Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. *Wien Med Wochenschr* (in German). 1997; 147(19-20): 463-4.
- Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus. *Microbiol Immunol.* 1994; 38(12): 951-5.

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
14. Referências

- Halstead SB. Selective primary health care: strategies for control of disease in the developing world: XI. Dengue. *Rev. Infect. Dis.* 1984; 6:251-264.
- Halstead SB. Pathogenesis of dengue: challenges to molecular biology. *Science* 1988; 239:476-481.
- Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). *Japanese J. Trop. Med. Hygiene* 1994; 22: 9-12.
- Lam SK. Dengue haemorrhagic fever. *Rev. Med. Micro.* 1995; 6:39-48.
- Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
- Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. *Clin Diagn Lab Immunol.* 2003 Jul; 10(4): 725-8.

- Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. *Wien Med Wochenschr (in German).* 1997; 147(19-20): 463-4.
- Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus. *Microbiol Immunol.* 1994; 38(12): 951-5.

Rev. 0, 2020-10-19 JC

Symbol	Deutsch	English	Français	Español	Italiano	Polski
	CE Konformitätszeichen	CE marking of conformity	Conforme aux normes européennes	Conformidad europea	Conformità europea	Znak zgodności CE
	Gebrauchsanweisung beachten	Consult instructions for use	Consulter la notice d'utilisation	Consúltense las instrucciones de uso	Consultare le istruzioni per l'uso	Przestrzegać instrukcji obsługi
	<i>in-vitro</i> -Diagnostika	<i>in-vitro</i> diagnostic medical device	Dispositif médical de diagnostic <i>in-vitro</i>	Producto sanitario para diagnóstico <i>in-vitro</i>	Dispositivo medico- diagnostico <i>in-vitro</i>	Tylko do diagnostyki <i>in-vitro</i>
	Temperaturbegrenzung	Temperature limitation	Limites de température	Límite de temperatura	Limiti di temperatura	Temperatura przechowywania
	Chargenbezeichnung	Batch code	Numéro de lot	Código de lote	Codice lotto	Numer serii
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
	Verwendbar bis	Use by	Utiliser jusqu'au	Fecha de caducidad	Utilizzare entro	Data ważności
	Bestellnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Riferimento di Catalogo	Numer katalogowy
	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
	Ausreichend für <n> Ansätze	Sufficient for <n> tests	Suffisant pour "n" tests	Suficiente para <n> utilizaciones	Sufficiente per "n" saggi	Wystarczający na <n> Powtórzeń

Symbol	Português	Český	Suomi	Svenskt	Nederlands	Dansk	Norsk
	Conformidade com as normas europeias	CE certifikát	CE-merkitty	CE-märkning	CE-markering	CE-mærkning	CE standardisert
	Consultar as instruções de utilização	Viz návod k použití	Katso käyttöohjetta	Läs bruksanvisningen	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen	Les bruksanvisning nøye
	Dispositivo médico para diagnóstico <i>in-vitro</i>	Diagnostický zdravotnický prostředek <i>in-vitro</i>	<i>in-vitro</i> - diagnostiikkaan tarkoitettu lääkinnällinen laite	Medicinteknisk produkt avsedd för <i>in-vitro</i> -diagnostik	Medisch hulpmiddel voor <i>in-vitro</i> diagnostiek	Medicinsk udstyr til <i>in-vitro</i> -diagnostik	<i>in-vitro</i> diagnostic medisinsk enhet
	Limites de temperatura	Teplotní omezení	Lämpötilarajat	Temperatur- begränsning	Temperatuurlimiet	Temperatur- begrænsning	Temperatur begrænsning
	Código do lote	Kód šarže	Eräkoodi	Satsnummer	Code van de partij	Batchkode	Merkning
	Não reutilizar	Pro jednorázové použití	Kertakäyttöinen	Får inte återanvändas	Niet opnieuw gebruiken	Må ikke genbruges	Må ikke brukes om igjen
	Prazo de validade	Spotřebuje do	Käytettävä viimeistään	Används före	Houdbaar tot	Udløbsdato	Tidtaking
	Número de catálogo	Katalogov číslo	Luettelonumero	Listnummer	Catalogus nummer	Best il lingsnummer	Katalog nummer
	Fabricante	Výrobce	Valmistaja	Tillverkare	Fabrikant	Fabrikant	Produsent
	Suficiente para <n> test	Dostačuje pro <n> testů	Lukumäärä <n> test	Räcker till <n> test	Voldoende voor <n> test	Tilstrækkeligt til <n> test	Tilstrækkelig for<n> tester

Our Teams

Germany:

Regensburg

Tel: +49 941 290 10-0
Fax: +49 941 290 10-50

Moers

Tel: +49 2841 99820-0
Fax: +49 2841 99820-1

Austria:

Tel: +49 941 290 10-29
Free Tel: 0800 291 565
Fax: +49 290 10-50
Free Fax: 0800 298 197

UK & Ireland:

Tel: +49 941 290 10-18
Free Tel –UK: 0808 234 1237
Free Tel – IRE: 1800 555 080
Fax: +49 290 10-50

France:

France Tel: 0800 915 240
France Fax: 0800 909 493

Switzerland

Swiss Tel: 0800 564 720
Swiss Fax: 0800 837 476

Belgium

Belgium Tel: 0800 718 82
Belgium Fax: 0800 747 07

Luxembourg

Lux, Tel: 800 211 16
Lux, Fax: 800 261 79

Spain:

Tel: +49 941 290 10-759
Free Tel: 900 938 315
Fax: +49 941 290 10-50
Free Fax: 900 984 992

Italy:

Tel: +49 941 290 10-34
Fax: +49 941 290 10-50

Poland:

Tel: +49 941 290 10-44
Free Tel: 00 800 491 15 95
Fax: +49 941 290 10-50
Free Fax: 00 800 491 15 94

Portugal:

Tel: +49 941 290 10-735
Tel, Verde: 800 849 230
Fax: +49 941 290 10-50
Fax Verde: 800 849 229

Netherlands:

Tel: +31 30 75 600
Free Tel: 0800 0222 890
Fax: +31 70 30 30 775
Free Fax: 0800 024 9519

Nordic countries:

Denmark

Tel: +31 703075 605
Free Tel: 808 887 53

Finland

Tel: +31 703075 606
Free Tel: 0800 918 263
Free Fax: 0800 918 262

Norway

Tel: +31 703075 605
Free Tel: 800 16 731

Sweden

Tel: +31 703075 605
Free Tel: 020 79 09 06