

PLASTOD S.p.A.	Title	PLASTER ON SPOOL IN WHITE FABRIC, NON PINKED EDGES, WITH ZINC OXIDE ADHESIVE, LATEX FREE	Ed. 01 of 25/03/04	TECHNICAL DATA SHEET		
<i>QUALITY ASSURANCE</i>	Sect.	DESIGN AND DEVELOPMENT	Rev. 1 of 06/04/06	06	005	0

Backing

Composition	100% cotton
Colour	white
Weight (g/m²)	90-110
MD Tensile Strength (N/cm)	> 40

References

Supplier's data sheet

IST 26 Sez. 7.5

IST 23 Sez. 7.5

Adhesive

Composition	white zinc-oxide artificial thermoplastic rubber, latex free
Weight (g/m²)	55-65

Supplier's data sheet

IST 26 Sez. 7.5

Coated backing

Weight (g/m²)	145-175
Adhesion to steel by peeling a 180° (N/cm)	> 2,5
M.V.T.R. (g/m²/24h/37°C)	< 500 (non perforated)
Water resistant	no

IST 26 Sez. 7.5

IST 21 Sez. 7.5

IST 27 Sez. 7.5

Finished product

Core	plastic spool
Edges	non pinked

Tolerances on final dimensions

Length (%):

≤ 3 m	± 2
> 3 m	± 2

Width (mm):

≤ 50 mm	± 1,5
> 50 mm	± 2

Joints per plaster

≤ 3 m	none
> 3 m	1

*Maximum percentage of
plasters with joints*

10

Legenda:

N = Newton = 0,1Kg
CD = cross direction

MD = machine direction

M.V.T.R.= Moisture Vapour Transmission Rate

DOCUMENT CONTROL		
<i>ISSUED BY: FRANCESCA MASSARIA</i>	<i>SIGN:</i>	<i>DATE: 06/04/06</i>
<i>VERIFIED BY: GIULIANO MICHELINI</i>	<i>SIGN:</i>	<i>DATE: 06/04/06</i>
<i>APPROVED BY: FRANCESCA MASSARIA</i>	<i>SIGN:</i>	<i>DATE: 06/04/06</i>

signed in original