

Cederroth Eye & Wound Cleansing spray, ref. 100EFA100



General product description

- A portable spray bottle filled with 150 mL physiological (0.9%) unbuffered sodium chloride solution, preservative free and sterile. The nozzle is suitable for use on both eyes and wounds. The actuator is protected by a plastic cap and a tamper seal. The nitrogen propellant is not in contact with the solution. The solution remains sterile until last drop. The Cederroth Eye & Wound Cleansing spray is intended to be used for rinsing and cleansing out dust and dirt from eyes and wounds.
- Storage information : store in temperatures between 0°C and +40°C.

Detailed product information

- Content per consumer unit
- **Material:** - Liquid: Sterile 0.9% Sodium Chloride (NaCl). Preservative free ;
 - Can: Aluminum ;
 - Bag-on-valve: 4 layer PP-OPA-ALU-PET ;
 - Actuator: White POM ;
 - Cover CAP: Transparent PP.
- **Shelf life:** 5 years from production date

Logistic information

- GTIN: 7310617260005.
- Content: 150 ml.
- Size: W 53 x L 53 x H 162 mm.
- Gross weight: 191 g.