

Cederroth Burn Gel Dressing, ref. 100EGA100



General product description

- A foam dressing impregnated with burn gel for the immediate first aid treatment of minor burns, scalds and sunburns. The gel rapidly cools and soothes the injury, prevents further tissue damage and offers a rapid and ongoing pain relief.
- The dressings are individually packed in a sealed, terminally sterilized pouches, that easily can be teared open. Two Cederroth Burn Gel Dressings are included in each consumer box.
- The Cederroth Burn Gel Dressing box fits into Cederroth First Aid Stations.
- Storage information : dry storage in room temperature. For maximum protection, keep in original package until use.

Detailed product information

- Content per consumer unit.
- **Product/Material:** Burn gel dressing in sealed pouch, sterile. Single use.
- **Dimensions:** - Dressing: 100 x 100 mm, quantity: 2 ;
- Pouch: 124 x 124 mm, quantity: 1.
- **Material:** - Burn Gel: Water(96%) based gel with propylene glycol, glycerol, emulsifier and surfactant ;
- Dressing: Polyurethane foam, 3 mm thick ;
- Pouch: Laminate of PE/ALU/PET ;
- Box: Carton.
- **Shelf life:** 5 years.

Logistic information

	CU	SKU	Pallet
GTIN	7310619019007	7310809019008	7312029019004
Content	2 dressings	24 CU	14 SKU/layer x 7 layers
Size (w x l x h mm)	126 x 36 x 65	260 x 205 x 180	Pallet 800 x 1200 x 1410
Gross weight	95 g	2433 g	263 kg
Packing material	Board 13 g	Corrugated board 153g	---